



CLEANING, DISINFECTING AND STORAGE SOLUTION

Please read carefully and keep this package insert for future use in case you have a problem.

DIRECTIONS FOR USE

GENERAL INSTRUCTIONS:

- Always wash and rinse your hands before handling your contact lenses.
- Clean and rinse the right lens first (always the same lens to avoid mix-ups).
- Optimum CDS is **NOT FOR USE DIRECTLY IN THE EYE.**

CLEAN:

Upon removing the right lens from the eye, place the lens in the palm of the hand and cover the lens completely with 10 drops of Optimum CDS. Rub the lens for 30 seconds with the index finger of the other hand creating a thick sudsy pool of solution.

NOTE: The greater the sudsing action the better the cleaning performance.

NOTE: An enzyme cleaner may also be used if recommended by your eye care practitioner.

RINSE:

Hold the lens under fresh running tap water. Be sure the drain of the sink is closed.

DISINFECT/STORE:

Place the lens in the right chamber of your lens storage case and completely cover the lens with Optimum CDS. Close the lid tightly and store the lens in the unopened storage case for at least six (6) hours. Repeat the above procedure for the left lens.

HOW SUPPLIED:

Optimum CDS is supplied sterile in 4 fl. oz. (118mL) plastic bottle. Bottles are marked with lot number and expiration date.

DESCRIPTION:

Optimum CDS is a sterile solution containing the lauryl sulfate salt of imidazoline, octylphenoxypolyethoxyethanol and preserved with benzyl alcohol 0.3% and disodium edetate 0.5%.

ACTIONS:

Optimum CDS is a combination cleaning, disinfecting and storage solution. Optimum CDS loosens and removes accumulations of film, debris and surface deposits from the lenses, and destroys harmful microorganisms on the surface of fluorosilicone acrylate and silicone acrylate GP contact lenses and hard contact lenses.

INDICATIONS:

Optimum CDS is indicated for use in the cleaning, chemical disinfection, and storage of fluorosilicone acrylate and silicone acrylate GP contact lenses and hard contact lenses.

CONTRAINDICATIONS:

Please note carefully the ingredients as listed. If you are allergic to any ingredient in Optimum CDS, do not use this product

WARNINGS:

- To avoid contamination, do not touch dropper tip to any surface.
- Close cap tightly after each use.
- Optimum CDS is **NOT FOR USE DIRECTLY IN THE EYE.**
- Optimum CDS must be completely rinsed from the lenses after disinfection and prior to wetting and insertion.
- If accidentally instilled in the eye, the eye will sting and become red. Flush the eye repeatedly with cool tap water. If pain or redness persists, consult a physician. If swallowed, drink extra water.
- Lens care procedures as recommended by your eye care practitioner must be followed. Failure to follow these procedures may result in serious eye infections. If any unexplained eye discomfort, watering, vision change, or redness of the eye occurs, immediately remove your lenses and consult your eye care practitioner to identify the cause and begin necessary treatment.

PRECAUTIONS:

- Always wash and rinse your hands before handling your lenses.
- Always use fresh Optimum CDS in your storage case.
NOTE: After inserting your lenses, always empty your lens storage case, rinse, and allow to dry.
- Never re-use the solution in your lens storage case.
- Store at room temperature 15-30°C (59-86°F).
- Keep out of reach of children
- Use before the expiration date stamped on carton and bottle label.
- Never heat Optimum CDS Solution on your GP lenses.
- Optimum CDS Solution is **NOT FOR USE WITH SOFT (HYDROPHILIC) CONTACT LENSES.**

ADVERSE REACTIONS

(PROBLEMS AND WHAT TO DO):

The following problems may rarely occur while wearing contact lenses:

- Redness of the eye.
- Eyes stinging, burning, or itching.
- Excessive watering (tearing) of the eyes.
- Unusual eye secretions.
- Reduced sharpness of vision (visual acuity).
- Blurred vision.
- Sensitivity to light (photophobia).
- Dry eyes.

If you notice any of the above problems, IMMEDIATELY remove and examine your lenses. If the problem stops and the lenses appear to be undamaged, thoroughly clean, rinse, and disinfect the lenses and reinsert them. If the problem continues, or a lens appears to be damaged, IMMEDIATELY remove your lenses and consult your eye care practitioner. DO NOT reinsert a damaged lens. If any of the above problems occur, a serious condition such as infection, corneal ulcer, neovascularization, or iritis may be present. Seek immediate professional identification and treatment of the problem to avoid serious eye damage. See your *Instructions for Wearers* Booklet for more information.

**EACH CONTAINER IS TAMPER-EVIDENT SEALED.
IF SEAL AROUND THE BOTTLE IS MISSING OR
BROKEN, DO NOT USE THIS PRODUCT.**

Distributed by:
Lobob Laboratories, Inc.
San Jose, CA 95131
www.LobobLabs.com

Comments or questions:
1(800) 835-6262
Info@LobobLabs.com